EN_PERSONAL BODY SCALE

THIS DEVICE IS DESIGNED TO MONITOR YOUR WEIGHT OVER TIME, ENSURING YEARS OF RELIABLE SERVICE. THE MANUAL PROVIDES DETAILED INSTRUCTIONS FOR AN OPTIMAL EXPERIENCE. THANK YOU FOR CHOOSING OUR SCALE!

WARNINGS

Read the instructions before use to avoid any danger due to misuse and keep it for further use in the future.

- Place your scale on a hard flat surface will ensure the greatest accuracy.
- Always weigh yourself on the same scale placed on the same floor surface. Uneven floors may affect the readings.
- Clean your scale with a damp cloth but do not allow water to get inside the scale, this can damage the electronics components.
- Do not use any chemical clearing agents.
- Do not store the scale in an upright position when not in use as this will drain the batteries.
- Treat your scale with care as it is a precision instrument, do not drop it or jump up or down on it.
- This product is intended only for the use for which it was designed, namely as an electronic scale for household use only. Any other use is to be considered improper and therefore dangerous.
- This scale is not to be used to determine the weight of objects or substances in commercial transactions, for the manufacture of medicines, for calculating tolls, rates, fees, taxes, premiums, fines, remuneration, allowances or fees of a similar type determined by weight.
- This scale may be used to keep track of body weight but should not be used to make diagnoses or define medical treatments. Values that may seem unusual or out of the norm should always be discussed with your physician. In order to avoid danger to the user, you should never change the dosage of any medication or therapy prescribed by your doctor.
- The product is not a toy.
- Careful supervision necessary in the presence of children, persons with disabilities or those not trained in its use.
- Product intended for indoor use. Do not use outdoors or expose to the weather.
- Use in dry and dust-free places.
- Do not use the device near bathtubs, showers, swimming pools.
- Never immerse the product in water or other liquids.
- Keep the product away from heat sources, hot and/or flammable surfaces.
- Check the condition of the product periodically and contact your dealer or qualified personnel if there are any abnormalities.
- Never disassemble the product or attempt to open it as this will void the warranty.
- When the scale is not in use, do not leave weights resting on the platform.
- Remove the batteries if the scales are not going to be used for a long period of time.
- Batteries cannot be disposed of with normal household waste. Find out about local regulations regarding the separate collection of batteries: proper disposal avoids negative environmental and health consequences.

SPECIFICATIONS

- Personal scale (LB/KG/ST)
- 6mm tempered glass platform
- Capacity: 180kg/396lb/28st 4lb
- Graduation: 100g (0.2lb/1/4st)
- Size: 300x300x21mm
- LED display with white digital; SEGMENT HEIGHT 28mm
- LED size: 33x78mm
- Equipped with high precision strain gauge sensor
- Step on & auto off
- Overload and low battery indication
- Power: 3x1.5V AAA battery (NOT INCLUDED)

POWER SUPPLY

- Power: 3xAAA Battery
- Please remove the insulation sheet from the battery contacts before using. Replace the batteries when the low battery indicator shows.

OPERATING INSTRUCTIONS

- **STEP 1.** Press the platform centre and remove your foot.
 - **STEP 2.** "0.0" will be displayed
 - **STEP 3.** The scale will switch off and now it is ready for use.
- This initialisation process must be repeated if the scale is moved. At all other times step straight on the scale.
- **STEP 4.** Gently step on the scale. Put your feet evenly on the scale without moving and wait whilst the scale compute weight. Your correct weight will then be displayed in the following way: 80.3 Kg.

AUTOMATIC SWITCH OFF

When you step off the scale, it will automatically switch off. (If no weighting is carried out while the zero is being displayed, the scale will automatically switch off after approximately 6 seconds).

OVERLOAD INDICATOR

When the scale is over the maximum capacity, the following will appear on the screen: ERR

BATTERIES

When the batteries need replacing the following will appear on the screen: LO

Please open the back board of the scale, then take the batteries carefully out of the battery compartment and insert the new ones following the correct polarity. Dispose of the old batteries in a safe place! Do not place on fire.

INFORMATION FOR USERS

pursuant to Legislative Decree N° 49 of 14 March 2014
"Implementation of Directive 2012/19 / EU on waste electrical and electronic equipment (WEEE)"
The crossed out bin symbol shown on the equipment indicates that the product must be at the end of its useful life collected separately from other waste. The user must, therefore, give the equipment complete with essential components reached the end of its life in suitable separate collection centers for electronic and electrotechnical waste, or return it to the retailer at the time of purchase of new equipment of an equivalent type, in the ratio of one to one, or 1 to zero for the equipment with longer side less than 25cm.
Adequate separate collection for the subsequent start-up of the discarded equipment for recycling, treatment and disposal compatible environment helps to avoid possible negative effects on the environment and health and promotes the recycling of materials which the equipment is composed of. Illegal disposal of the product by the user involves the application of sanctions administrative provisions pursuant to Legislative Decree no. Legislative Decree N°49 of 14 March 2014.

©2025 AAAmaze®. Produced by SEECON GROUP LIMITED - 88 Kwai Cheong Road, Unit 1203, 12/F., NO. 88 Kwai Cheong Road, Kwai Chung Hong Kong - info@seecongroup.com - Imported and distributed by HUBITECH Srl - Via Gian Galeazzo, 16 - 20136 Milano (Italy) on behalf of DIMO Spa
info@hubitech.it - info@aaamaze.com - All trademarks belong to their respective owners.

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©2025 AAAMaze®. Produced by SEECON GROUP LIMITED - 88 Kwai Cheong Road, Unit 1203, 12/F., NO. 88 Kwai Cheong Road, Kwai Chung Hong Kong - info@seecongroup.com - Imported and distributed by HUBITECH Srl - Via Gian Galeazzo, 16 - 20136 Milano (Italy) on behalf of DIMO Spa info@hubitech.it - info@aaamaze.com - All trademarks belong to their respective owners.